



THE UNIVERSITY *of* EDINBURGH

Edinburgh Research Explorer

Comparison of Doppler ultrasonic and oscillometric devices (with or without 2 proprietary optimizations) for non-invasive blood pressure measurement in conscious 3 cats

Citation for published version:

Cerna, P, Archontakis, PE, Cheuk, HOK & Gunn-Moore, D 2020, 'Comparison of Doppler ultrasonic and oscillometric devices (with or without 2 proprietary optimizations) for non-invasive blood pressure measurement in conscious 3 cats' *Journal of Feline Medicine and Surgery*.
<https://doi.org/10.1177/1098612X20932407>

Digital Object Identifier (DOI):

[10.1177/1098612X20932407](https://doi.org/10.1177/1098612X20932407)

Link:

[Link to publication record in Edinburgh Research Explorer](#)

Document Version:

Peer reviewed version

Published In:

Journal of Feline Medicine and Surgery

General rights

Copyright for the publications made accessible via the Edinburgh Research Explorer is retained by the author(s) and / or other copyright owners and it is a condition of accessing these publications that users recognise and abide by the legal requirements associated with these rights.

Take down policy

The University of Edinburgh has made every reasonable effort to ensure that Edinburgh Research Explorer content complies with UK legislation. If you believe that the public display of this file breaches copyright please contact openaccess@ed.ac.uk providing details, and we will remove access to the work immediately and investigate your claim.



Comparison of Doppler ultrasonic and oscillometric devices (with or without proprietary optimizations) for non-invasive blood pressure measurement in conscious cats

Petra Cerna^{1,2}, Panos E Archontakis³, Hester O K Cheuk⁴, Daniëlle A Gunn-Moore^{1,5}

¹ The Royal (Dick) School of Veterinary Studies, University of Edinburgh, UK

² The University of Veterinary and Pharmaceutical Sciences Brno, CZ

³ Department of Clinical Science and Services, Royal Veterinary College, University of London, UK

⁴ First Vets Ltd, Newcastle Upon Tyne, UK

⁵ The Roslin Institute, University of Edinburgh, UK

Corresponding author: MVDr. Petra Cerna MRCVS petra.cerna.mvdr@gmail.com

MVDr. Petra Cerna MRCVS, Department of Internal Medicine, Hospital for Small Animals, The Royal (Dick) School of Veterinary Studies, University of Edinburgh, UK

Panos E Archontakis DVM, MRCVS, Department of Clinical Science and Services, Royal Veterinary College, University of London, UK

Hester O K Cheuk BVM&S MRCVS, First Vets Ltd, Newcastle Upon Tyne, UK

Prof. Daniëlle A Gunn-Moore BSc(Hon), BVM&S, PhD, MANZCVS, FHEA, FRSB, FRCVS, RCVS Specialist in Feline Medicine, Department of Internal Medicine, Hospital for Small Animals, The Royal (Dick) School of Veterinary Studies, University of Edinburgh, UK

Abstract:

Objectives: This study compared Doppler and oscillometric (PetMAP+) devices (with or without proprietary optimizations) for the non-invasive measurement of blood pressure in conscious cats.

Methods: Twenty-three cats were enrolled; however, five were excluded as fewer than five measurements were obtained for each assessment. All measurements were obtained according to American College of Veterinary Internal Medicine (ACVIM) consensus guidelines. Oscillometric device modes A and B were operated according to manufacturer's guidelines. Doppler and oscillometric devices were alternatingly used first.

Results: Systolic arterial blood pressure (SAP) measurements were obtained by Doppler (SAPd) and oscillometry; the mean of each set of five values was used for statistical analysis. There was a significant difference between SAPd and SAP measurements on oscillometric mode A ($p<0.001$) and B ($p<0.001$). While both modes measured SAP higher than SAPd, B had a smaller bias (+15.72mmHg) and narrower limits of agreement. There was also a significant difference between SAPd and mean arterial pressure (MAP) on oscillometric modes A ($p=0.002$) and B ($p<0.001$). Both modes' MAP readings were lower than SAPd, oscillometric

A MAP was closer to SAPd (-14.94mmHg), with a smaller bias and narrower limits of agreement.

Conclusions and relevance: The findings support that Doppler and oscillometric devices cannot be used interchangeably, with or without proprietary optimizations. Methodology should always be taken into account, and reference intervals need to be defined for the different methodologies. Until methodology-specific reference intervals are published, definitive diagnosis of hypertension and sub-staging of patients with kidney disease according to the IRIS guidelines remains challenging.

Introduction:

Blood pressure (BP) is a valuable measurement in feline veterinary medicine. Non-invasive BP measurement in conscious cats is a vital part of routine monitoring of older individuals, and for patients on medication that may affect blood pressure. Moreover, it is an important indicator of cardiovascular function in critically ill¹ and anaesthetised animals.² Monitoring cases of known systemic hypertension and assessing cats with potential causes of secondary hypertension are major indications for routine BP measurements in middle aged to older cats (from nine years of age, onwards).³⁻⁵ Primary, also called idiopathic, hypertension has been reported in 13-20% of cats with hypertension.^{4,6,7} It is not yet fully understood if some of these cats might have non-azotaemic chronic kidney disease (CKD) or if environmental factors may play a role. Genetic predisposition to primary hypertension has been reported in humans⁸; however, this has not yet been documented in cats. Secondary hypertension is strongly associated with diseases such as CKD, hyperthyroidism, primary hyperaldosteronism and, less commonly, other rare diseases such as hyperadrenocorticism and pheochromocytoma.^{4,9,10} The most common condition associated with feline hypertension in cats is CKD. Azotaemia and hypertension are intimately linked; azotaemia has been found in up to 74% of hypertensive cats, and 19-65% of cats with CKD have been found to be hypertensive.⁹⁻¹² However, the prevalence and severity of hypertension does not appear to be related to the severity of the CKD.^{9,13,14} Hypertension is most likely to cause disease in tissues with a rich arteriolar supply^{3,15} and an early recognition of hypertension is required to prevent irreversible organ damage to the heart, brain, eyes, and kidneys, as these are particularly vulnerable to hypertensive injury or so called target organ damage (TOD).^{3,4,7,9,10,16}

Blood pressure can be measured directly via intra-arterial catheterization or indirectly via non-invasive devices that incorporate a compressive cuff placed on a limb or the base of the tail. Invasive blood pressure (IBP) measurement is recognized as the gold standard; however, this method is not commonly used as it is impractical in routine clinical practice; it requires technical expertise, advanced medical equipment and is uncomfortable for the patient.^{3,17} For this reason, indirect BP measurements obtained via non-invasive blood pressure (NIBP) devices such as Doppler ultrasonic sphygmomanometry or oscillometry are more commonly used.

Doppler ultrasonic and oscillometric machines are popular NIBP devices currently available in the UK. Doppler devices have been extensively used in feline medicine,^{18,19} with one study demonstrating good correlation and accuracy compared with direct BP assessment.²⁰ They detect blood flow utilizing the Doppler Effect on moving red blood cells. Systolic (SAP) and diastolic (DAP) arterial BP can be measured; however, the DAP measurement can be difficult to obtain, and the accuracy of these measurements is dependent on the operator's experience, particularly when measuring DAP.^{2,21} It can be difficult to obtain BP values in fractious and uncooperative cats. Oscillometric devices are technically easier to use than Doppler machines. They rely on the detection of oscillations in the artery and built-in algorithmic calculations to produce automated SAP, mean arterial pressure (MAP) and diastolic arterial pressure (DAP) measurements. It is important to remember that while Doppler devices all use the same methodology, each oscillometric device uses the particular algorithm devised by its manufacturer. In addition, while they are usually easier to use than Doppler devices, it is still not always possible to achieve BP measurements using them.^{2,20,22-26}

Studies have tried to establish reliable reference intervals for SAP in healthy conscious and anaesthetised cats using intra-arterial, Doppler and oscillometric equipment (Table 1); a number of studies have suggested that traditional oscillometry measurements are less accurate than Doppler in conscious cats, and often underestimate SAP at higher values.^{2,20,22,24-26} One study found no significant difference between mean SAP readings when they compared a Doppler device (CAT Doppler, Thames Medical, UK) and an oscillometric device (Memoprint SCB, medVetGmbH, Germany); however, the oscillometric machine produced significantly higher estimates of DAP.²³ When looking at Table 1, the mean SAP by IBP ranged from 118 to 160 mmHg, compared to 115 to 139 mmHg for oscillometric devices, and 121 to 162 mmHg for Doppler devices, so there is a great deal of variation, irrespective of the methodology or device used.

It is important to remember that, as in other species, BP in cats varies markedly within and between individuals.^{19,27-29} The cat's activity level, excitement and anxiety can play a significant role in the accurate measurement of BP. It is important to give cats time to calm down and acclimatise after entering the consultation room and to use as stress-free approach as possible.²⁷⁻²⁹ One study showed up to an 80 mmHg change in SAP in response to a simulated clinic visit, showing the potential magnitude of situational hypertension in healthy cats.²⁷ Systolic arterial BP can be affected by many external variables including the operator, conditions within the environment when and where the BP measurement is taken, the equipment, the position of the cat, the size of the cuff, and the site of measurement.^{19-21,23,27,29-33} Accepting the difficulties of gaining reliable BP readings, once each clinic establishes its own reliable protocol and its own bias on the reference interval (e.g. being higher in a noisy clinic), BP can be reasonably reliably assessed in pet cats.

While there appear to be no gender or breed effects on BP,^{34,35} a large longitudinal study established a small but significant increase in BP with increasing age, equating to ~1-2 mmHg per annum for cats over nine years of age.¹⁴ Indirect SAP measurements exceeding 160 mmHg have been associated with TOD in several studies.^{4,7,10,12,36-39}

The purpose of this prospective study was to evaluate the agreement of SAP measured using a Doppler machine (Parks model 811, Parks Medical Electronics, Inc., Oregon, USA) and an oscillometric device (PetMAP+, model 7142, Ramsey Medical, Inc., Florida, USA) in two modes: mode A (PetMAP+ mode for proprietary optimisations [PPO]) and mode B (Optimized None [ON]). The PPO mode has built-in proprietary algorithms designed to optimise the estimation to equate to intra-arterial pressure. As IBP measurements could not be performed as a comparison, the accuracy of either machine could not be compared against a gold standard control. The hypothesis was that there would be good agreement between SAP measurement by the Doppler and the oscillometric device.

Methods:

Animals:

The study was conducted in the Hospital for Small Animals, of the Royal (Dick) School of Veterinary Studies, in 2015 after the veterinary ethical review committee approval of this study. Feline in-patients and healthy cats belonging to staff members were recruited for this study. Cats that were fractious, uncooperative and moving excessively were excluded. Any cat

receiving a drug known to have an effect on BP, or that had been anaesthetised or sedated within the previous 12 hours was also excluded.

Blood Pressure measurements:

To achieve standardized conditions all BP measurements were obtained according to the guidelines in the American College of Veterinary Internal Medicine (ACVIM) consensus statement on feline hypertension.³ A single well-trained operator with the help of a consistent assistant made all BP measurements in a designated Cat Ward. The Doppler (Parks model 811, Parks Medical Electronics, Inc., Oregon, USA) and oscillometric (PetMAP+, model 7142, Ramsey Medical, Inc., Florida, USA) devices were used alternatively as the first device. Both devices were operated carefully according to the manufacturer's guidelines. Measurements were taken on the antebrachium of the non-weight-bearing forelimb, with the cat sitting comfortably and the limb positioned at the level of the heart. The same inflatable cuff provided by Ramsey Medical, Inc. was used for both devices. An ideal cuff width between 30-40% of limb circumference was chosen. Prior to all BP measurements, each cat was given a 5-minute acclimatization period to the inflatable cuff in the aim of reducing stress-induced situational hypertension. All BP measurements were performed successively with the first measurement of each device discarded and the subsequent five measurements averaged.

Doppler measurement: the hair over the probe site was not clipped to reduce stress to the cat. Coupling gel was placed over the palmar aspect of the common digital artery and the Doppler probe was applied to obtain a clear audible sound. The cuff was manually inflated until the pulse signal was no longer audible, then gradually deflated. The reappearance of flow sounds in conjunction with the manometer reading represented the SAP (since it was for the Doppler device; SAPd). The cuff was completely deflated between each measurement.

The oscillometric measurement: the device was programmed to Clinic Mode for full operator control over the timing and frequency of BP measurements. Measurements were taken from the ON mode and PPO mode for "Cat Forearm" to obtain SAPon and SAPppo. Indirect MAP and heart rate readings were taken automatically, from which SAP and DAP were calculated by the device. Readings that were obtained during a sudden movement or displayed a poor oscillometric envelope were discarded and subsequently replaced.

Statistic Method:

The mean of the five SAP measurements was used for statistical analysis. Bland Altman analysis was used to determine the agreement between SAP measurements of the Doppler and SAP of the oscillometric device in both modes, as well as SAP measurements of the Doppler and MAP of the oscillometric device in both modes. The differences in SAP (SAPon – SAPd) or (SAPppo – SAPd) were plotted against the mean of each pair of measurement. The bias was defined by the mean differences between the two methods, (Σ [SAPon – SAPd]/2) and (Σ [SAPppo – SAPd]/2). The precision was calculated from the standard deviation (SD) of the differences and limits of agreement were calculated from the bias \pm 1.96 x SD. A paired sample T-test was run to determine whether there was a statistically significant mean difference in the SAPd compared to the SAP measured by the two oscillometric device modes. Results were presented as means \pm SD. A P-value less than 0.05 was considered significant. Overall statistical analysis was achieved using GraphPad Prism version 8.2.1.

Results:

Successful Readings

Twenty-three cats were enrolled; however, five cats were excluded as fewer than five BP measurements could be obtained for each assessment. In total, 18 cats were included in the study; 10 males and eight females, with a median age of eight years (range nine months to 17 years).

A total of 270 measurements were taken from the 18 conscious cats, with 90 readings from the Doppler, and each of the oscillometric device modes (Table 2). Measurements for these 18 cats were successfully obtained in 100% of attempts using the Doppler device and 94.4% of attempts using the oscillometric device. Five of the oscillometric readings were unsuccessful as they displayed a poor oscillometric waveform envelope and were therefore discarded as unreliable results. These are not present in the data in Table 3. The hypertensive patients were followed up and treated appropriately.

Agreement between SAPd and SAP (oscillometric device)

Bland-Altman analyses of the differences between SAP measurements from the two devices were plotted against the mean of the SAPd (Figure 1 and 2). Comparison of oscillometric mode A to SAPd gave an average bias of +21.28 mmHg (95% CI; 12.49 to 30.07 mmHg). The lower limit of agreement (LLA) was calculated as -13.36 mmHg, and the upper limit of agreement (ULA) was calculated as +55.92 mmHg; the difference between the upper and lower limit was 69.28 mmHg. Comparison of oscillometric mode B with SAPd produced an average bias of

+15.72 mmHg (95% CI; 7.46 to 23.98 mmHg). The limits of agreements were calculated as -16.83 for the LLA and +48.28 for the ULA; the difference between LLA and ULA was 65.11 mmHg.

There was a significant difference between SAPd and oscillometric mode A ($p<0.001$) and B ($p<0.001$). As indicated by the bias, both oscillometric device functions overestimated the SAP compared to SAPd. However, the oscillometric mode B performed better for the detection of SAP as demonstrated by the smaller bias and narrower limits of agreement. A possible proportional bias was noticed on both modes as the degree of the over-estimation was greater at higher compared to lower BP measurements: this suggests there is less agreement between the methods when measuring higher BP.

Agreement between SAPd and MAP (oscillometric device)

The differences between MAP measurements were plotted against the SAPd (Figure 3 and 4). Comparison of MAP on oscillometric mode A with the SAPd gave an average bias of -14.94 mmHg (95% CI: -23.66 to -6.22 mmHg). The limits of agreements were calculated as -48.18 mmHg for the LLA and +18.30 mmHg for the ULA; the difference between LLA and ULA was 66.48 mmHg. The comparison of MAP on oscillometric mode B with the SAPd produced an average bias of -20.24 mmHg (95% CI: -26.78 to -13.70 mmHg). The limits of agreements were calculated as -45.17 for LLA and +4.70 for ULA; the difference between ULA and LLA was 49.86 mmHg.

There was a significant difference between the SAPd and oscillometric MAP in mode A ($p=0.002$) as well as between SAPd and oscillometric MAP in mode B ($p<0.001$), with the MAP being lower than SAPd. Oscillometric mode A was closer to the SAPd as demonstrated by the smaller bias and narrower limits of agreement. A possible proportional bias was noticed on mode A as the degree of the under-estimation was greater at lower than at higher BP measurements. This trend was not noticed with mode B.

Discussion:

The aim of this study was to compare the Doppler BP device (Parks model 811) to the more recently marketed veterinary oscillometric machine (petMAP+). The study evaluated the bias and precision of SAP measurements with this oscillometric device set in two different modes; petMAP+ PPO (mode A) and ON (mode B). The study showed that the oscillometric device

had poor agreement with paired Doppler measurements in conscious cats. The reported positive bias and wide variability of SAP for both oscillometric modes indicated poor agreement and precision. In both modes, the mean SAP was higher than that for the Doppler (bias of +21.28 mmHg and +15.72 mmHg for mode A and B, respectively). It has been suggested that SAPd may be more reflective of MAP in cats.⁴⁰ To investigate this the oscillometric device MAP measurement was compared to SAPd and found to be lower (bias of -14.94 mmHg and -20.24 mmHg for mode A and B respectively).

Other oscillometric devices are also reported to show significant variability in BP measurements; these included the Memoprint SCB (medVet GmbH, Germany),²³ Memo Diagnostic HDO Pro SpB (medVet GmbH, Germany)⁵ and Dinamap Model 8300 (Critikon, USA).²⁰ A study evaluating the previous petMAP device reported poor agreement between the indirect and direct BP methods, and concluded that the device could not be recommended for BP measurements in cats.⁴¹

The newer petMAP+ oscillometric device was evaluated in an unpublished study of 10 cats, and reported overall improvements to BP accuracy and oscillometric sensitivity; however, the cats were lightly anaesthetised.⁴² In the present study, using the same device as above, but in conscious cats, the readings with mode B (i.e. mode ON) showed a better accuracy and precision to SAPd than mode A (although it gave a mean bias of +15.72 mmHg). This is in agreement with the manufacturer, that states that petMAP+ ON should show better correlation with “other [indirect] BP devices”. However, the manufacturer states that PetMAP+ PPO (mode A) is designed to improve the correlation between NIBP to IBP measurements; it is also unique to the selected species (dog and cat) and cuff placement site (forearm, tail or hind-foot). Built-in proprietary algorithms have been formulated to produce readings that are 10-20% higher than SAPd, on the basis that Doppler devices consistently underestimate the intra-arterial SAP; hence it should reflect intra-arterial SAP more accurately. The tendency of Doppler devices to consistently underestimate intra-arterial SAP has been supported by several studies in anaesthetised cats.^{2,30,40} In the current study, using mode A, SAP gave a mean bias of +21.28 mmHg, so it was ~10-20% higher, *as per* the manufacturer’s claims. It has been previously described that SAPd is occasionally closer to MAP than SAP as measured by IBP.^{43,44} In the current study, MAP (by oscillometry) had a mean bias against SAPd of -14.94 mmHg in mode A, and -20.24 mmHg in mode B, which means MAP (by oscillometry) was overestimated by

the SAPd in both modes; however, the true relationship of MAP (by oscillometry) and MAP (by IBP) is unclear at this time.

Studies have assessed IBP, and compared it with oscillometric and Doppler devices, in conscious and anaesthetised cats (Table 1); unfortunately, no study has compared all three devices with both situations. Some studies have attempted to validate indirect BP devices in anesthetized patients, while fewer have been performed in conscious animals.^{19,20,22,23,34,47,48} In one study, the Doppler (Parks model 811, Parks Medical, USA) and oscillometric devices (Dinamap 8300, Critikon, USA) provided a reliable estimate of BP in anaesthetised cats; however, the Doppler device provided the most reliable measurements in conscious cats.²⁰ It has been suggested that traditional oscillometry is less accurate than Doppler in conscious cats as it often underestimates BP, especially at higher values and produces less reliable readings.^{2,20,23,24,26} In contrast, the current study found the oscillometric device overestimated the SAP compared to the Doppler. That said, its SAP on A mode was almost within the 10-20% the manufacturer defined, which may be closer to the real SAP; it also found MAP to be within the 10-20% the manufacturer defined in mode A. As with other studies, there was considerable variation between the two devices; although the mean of the Doppler device was close to those previously reported (from 121 to 162 mmHg; Table 1), the mean of the oscillometry device in both modes was higher than previously reported (from 115 to 139 mmHg; Table 1). Moreover, they were not defined enough to be usefully applied to ACVIM guidelines. Further studies are needed before devices that can meet the guidelines for measuring BP in conscious animals have been universally validated; ACVIM guidelines recommend that currently available devices should be used with a degree of caution.³

Consistent with several studies reporting measurement failures in the evaluation of oscillometric devices, albeit in anaesthetised cats,^{2,5,20} failure to obtain a reading occurred in approximately 5.6% of the measurement attempts with the oscillometric device in the current study. Previous user-specific errors that have been shown to have an adverse effect on BP determinations include the operator's experience with the Doppler device^{2,21} and improper cuff size in both Doppler and oscillometric devices.^{2,19} Since the same cuff was used for both machines, and the operator was experienced, these were unlikely to have played a significant role the current study's findings.

The main limitations of the study were its small sample size and that a comparison of the two non-invasive methods were not compared with a direct gold standard BP measurement. To compare both modes on the oscillometric device with the Doppler, each cat had to have over 15 individual BP assessments. This means that even some of the placid cats became restless, so a full assessment could not be performed, and cats had to be removed from the study. A larger number of cats would have improved the study, as would having more cats with blood pressure above and below the reference interval. However, since participation involved numerous BP measurements, the stress of inclusion would have been inappropriate for hypertensive patients. Most hypotensive cats are hospitalised in our Intensive Care Unit, where such a study is not allowed.

It has been suggested that lower reference intervals for normal SAP values should be used when using oscillometric devices, compared to those used with the Doppler, since two cats with ocular TOD were missed using the Memoprint oscillometric device (medVet GmbH, Germany).²³ However, the range of mean values found with both types of device is very variable (Table 1), and the current study found that oscillometric SAP was higher than SAPd, so this suggestion does not appear sensible. There is a real need for clarified reference intervals, one for SAPd and another for oscillometric devices, that is, if different oscillometric devices can be found to have the same reference interval (which has so far not proved to be the case; Table 1). Additional studies may be able to determine if it is possible to create a conversion calculation between Doppler and newer oscillometric devices. Given the potentially deleterious effects of missing hypertensive cats and over-treating cats falsely believed to be hypertensive, this area of feline medicine urgently needs more investigation.

Conclusions:

In conclusion, the oscillometric device did not produce SAP readings in good agreement with the Doppler device in conscious cats. While the petMAP+ PPO features did produce results 10-20% higher than SAPd, which the manufacture suggests is more consistent with IBP, there was still variability and a lack of gold standard. However, the beneficial qualities of the oscillometric device such as portability, intuitive user interface and minimal requirement of restraint still renders this device (and other oscillometric machines) a suitable choice in fractious and uncooperative cats. Due to the previously reported tendency of some oscillometric devices to underestimate SAP, especially at higher values,^{2,20,23-25} compared to the Doppler, while the current study found the PetMAP+ to overestimate SAP, all BP readings should be

integrated with the patient's history and physical examination. Importantly, so that changes over time can be determined, individual cats should be assessed using the same machine, cuff size and procedure each time (which should be recorded in their notes); this way, they become their own control, unrelated to any published reference intervals. Regardless of the BP measurement method and BP value, significant clinical signs that reflect TOD should be investigated via a fundoscopic examination, neurological evaluation and thoracic auscultation.¹⁶ Methodology-specific reference intervals are needed for more accurate sub-staging of patients with kidney disease and decision making about when to start antihypertensive treatment as the current variability makes decision making according to the IRIS, ISFM and ACVIM guidelines challenging.⁴⁹⁻⁵¹

Acknowledgements: The authors would like to thank all of the owners who allowed their cats to be used in this study, Dr Darren Shaw for his help with statistical analysis of this data and Luke Winter for cat holding assistance while the cats had their BP taken. The petMAP+ device was kindly lent by Ramsey Medical, Inc.

Conflict of Interest: The authors declared no potential conflicts of interest with respect to the research, authorship, and/or publication of this article.

Funding: The authors received no financial support for the research, authorship, and/or publication of this article.

Ethical Approval: The survey was approved by the Veterinary Ethical Review Committee (VERC) at the Royal (Dick) School of Veterinary Studies. The ethical review committee approved this study to be conducted under the legislation of Animals (Scientific Procedures) Act 1986.

Informed Consent: Informed consent (either verbal or written) was obtained from the owner or legal custodian of all animal(s) described in this work for the procedure(s) under taken. No animals or humans are identifiable within this publication, and therefore additional informed consent for publication was not required.

References:

1. Simpson KE, McCann TM, Bommer NX, et al. Retrospective analysis of selected predictors of mortality within a veterinary intensive care unit. *J Feline Med Surg* 2007; 9(5):364-8.
2. Binns SH, Sisson DD, Buoscio DA, et al. Doppler ultrasonographic, oscillometric sphygmomanometric, and photoplethysmographic techniques for noninvasive blood pressure measurement in anesthetized cats. *J Vet Intern Med* 1995; 9: 405–414.
3. Brown S, Atkins C, Bagley R, et al. Guidelines for the identification, evaluation, and management of systemic hypertension in dogs and cats. *J Vet Intern Med*. 2007; 21:542-558.
4. Elliot J, Barber PJ, Syme HM, et al. Feline hypertension: clinical findings and response to antihypertensive in 30 cases. *J Small Anim Pract* 2001; 42: 122-129.
5. Petric AD, Petra Z, Jerneja S, et al. Comparison of high definition oscillometric and Doppler ultrasonic devices for measuring blood pressure in anaesthetized cats. *J Feline Med Surg* 2010; 12:731-737.
6. Jepson RE, Elliott J, Brodbelt D, et al. Effect of control of systolic blood pressure on survival in cats with systemic hypertension. *J Vet Intern Med* 2007; 21: 402–409.
7. Maggio F, DeFrancesco TC, Clarke ET, et al. Ocular lesions associated with systemic hypertension in cats: 69 cases (1985 -1998). *J Am Vet Med Assoc* 2000; 217 (5): 695-70.
8. Oparil S, Zaman MA and Calhoun DA. Pathogenesis of hypertension. *Ann Intern Med* 2003; 139: 761–776
9. Kobayashi DL, Peterson ME, Graves TK, et al. Hypertension in Cats with Chronic Renal Failure or Hyperthyroidism. *J Vet Intern Med* 1990; 4(2): 58-62.
10. Littman MP. Spontaneous Systemic Hypertension in 24 Cats. *J Vet Intern Med* 1994; 8(2): 79-86.
11. Stiles J, Polzin D and Bistner SI. The prevalence of retinopathy in cats with systemic hypertension and chronic renal failure or hyperthyroidism. *J Am Anim Hosp Assoc* 1994; 30: 564–572.
12. Syme HM, Barber PJ, Markwell PJ, et al. Prevalence of systolic hypertension in cats with chronic renal failure at initial evaluation. *J Am Vet Med Assoc* 2002; 220: 1799–1804.
13. Karck J, von Spiessen L, Rohn K, et al. Interrelation between the degree of a chronic renal insufficiency and/or systemic hypertension and ocular changes in cats. [Article in German]. *Tierarztl Prax Ausg K Kleintiere Heimtiere* 2013; 41: 37–45.
14. Bijsmans ES, Jepson RE, Chang YM, et al. Changes in systolic blood pressure over time in healthy cats and cats with chronic kidney disease. *J Vet Intern Med* 2015; 29: 855–861.

15. Henik RA. Diagnosis and treatment of feline systemic hypertension. *Comp Cont Educ Pract* 1997; 19: 163–179.
16. Stephien RL. Feline systemic hypertension: diagnosis and management. *J Feline Med Surg* 2011; 13: 35-43.
17. Branson KR, Wasgner-Mann CC, Mann FA. Evaluation of an oscillometric blood pressure monitor on anaesthetized cats and the effect of cuff placement and fur on accuracy. *Vet Surg* 1997; 26: 347-353
18. Payne JR, Brodbelt DC, Luis FV. Blood pressure measurements in 780 apparently healthy cats. *J Vet Intern Med* 2017; 31: 15-21.
19. Sparkes AH, Caney SM, King MC, et al. Inter- and intra- individual variation in Doppler ultrasonic indirect blood pressure measurements in healthy cats. *J Vet Intern Med* 1999; 13: 314-318.
20. Haberman CE, Morgan JD, Kang CW, et al. Evaluation of Doppler Ultrasonic and Oscillometric Methods of Indirect Blood Pressure Measurements in Cats. *Int J Appl Res Vet M* 2004; 2(4): 279-289.
21. Gouni V, Tissier R, Misbach C, et al. Influence of the observer's level of experience on systolic and diastolic arterial blood pressure measurements using Doppler ultrasonography in healthy conscious cats. *J Feline Med Surg* 2015; 17(2): 94-100.
22. Martel E, Egner B, Brown SA, et al. Comparison of high- definition oscillometry – a non-invasive technology for arterial blood pressure measurement – with a direct invasive method using radio-telemetry in awake healthy cats. *J Feline Med Surg* 2013; 15: 1104–1113.
23. Jepson RE, Hartley V, Mendi M, et al. A comparison of CAT Doppler and oscillometric Memoprint machines for non-invasive blood pressure measurement in conscious cats. *J Feline Med Surg* 2005; 7: 147-152.
24. Sander C, Hörauf A and Reusch C. Indirect blood pressure measurement in cats with diabetes mellitus, chronic nephropathy and hypertrophic cardiomyopathy. *Tierarztl Prax Ausg K Kleintiere Heimtiere* 1998; 26: 110–118.
25. Möllenhoff A, Nolte I and Kramer S. Indirect blood pressure determination in cats using Doppler ultrasonic and oscillometric method. *Tierärztliche Praxis* 2001; 29: 191–197.
26. Pedersen KM, Butler MA, Ersbøll AK, et al. Evaluation of an oscillometric blood pressure monitor for use in anesthetized cats. *J Am Vet Med Assoc* 2002; 221: 646–650.
27. Belew AM, Barlett T and Brown SA. Evaluation of the whitecoat effect in cats. *J Vet Intern Med* 1999; 13: 134–142.

- 444 28. Mishina M, Watanabe N and Watanabe T. Diurnal variations of blood pressure in cats. *J*
445 *Vet Med Sci* 2006; 68: 243–248.
- 446 29. Slingerland LI, Robben JH, Schaafsma I, et al. Response of cats to familiar and unfamiliar
447 human contact using continuous direct arterial blood pressure measurement. *Res Vet Sci*
448 2008; 85: 575–582.
- 449 30. Grandy JL, Dunlop CI, Hodgson DS, et al. Evaluation of the Doppler ultrasonic method of
450 measuring systolic arterial blood pressure in cats. *Am J Vet Res* 1992; 53: 1166–1169.
- 451 31. Cannon MJ and Brett J. Comparison of how well conscious cats tolerate blood pressure
452 measurement from the radial and coccygeal arteries. *J Feline Med Surg* 2012; 14: 906–909.
- 453 32. Conti LMDC, Champion T, Guberman UC, et al. Comparison of indirect systolic blood
454 pressure on the forelimb and hindlimb of cats. *Rev Acad Ciênc Agrár Ambient* 2013; 11:
455 395–401.
- 456 33. Curtet JD, Busato A and Lombard CW. The use of memoprint in the cat. *Schweiz Arch*
457 *Tierheilkd* 2001; 143: 241–247.
- 458 34. Mishina M, Watanabe T, Fujii K, et al. Non-invasive blood pressure measurements in cats:
459 clinical significance of hypertension associated with chronic renal failure. *J Vet Med Sci*
460 1998; 60: 805–808.
- 461 35. Bodey AR and Sansom J. Epidemiological study of blood pressure in domestic cats. *J Small*
462 *Anim Pract* 1998; 39: 567–573.
- 463 36. Sansom J, Barnett KC, Dunn KA, et al. Ocular disease associated with hypertension in 16
464 cats. *J Small Anim Pract* 1994; 35: 604–611.
- 465 37. Bodey AR and Sansom J. Epidemiological study of blood between proteinuria, systemic
466 hypertension and glomerular filtration rate in dogs with renal and non-renal diseases. *Vet*
467 *Rec* 2008; 162: 141–147.
- 468 38. Wey AC and Atkins CE. Aortic dissection and congestive heart failure associated with
469 systemic hypertension in a cat. *J Vet Intern Med* 2000; 14: 208–213.
- 470 39. Kyles AE, Gregory CR, Wooldridge JD, et al. Management of hypertension controls
471 postoperative neurologic disorders after renal transplantation in cats. *Vet Surg* 1999; 28:
472 436–441.
- 473 40. Caulkett NA, Cantwell SL, Houston DM: A comparison of indirect blood pressure
474 monitoring techniques in the anesthetized cat, *Vet Surg* 27:370,1998.
- 475 41. Acierno MJ, Seaton D, Mitchell MA, et al. Agreement between directly measured blood
476 pressure and pressures obtained with three veterinary-specific oscillometric units in cats. *J*
477 *Am Vet Med Assoc* 2010; 237(4): 402–406.

42. Muir W. Accuracy of a species and cuff site optimised oscillometric BP monitor in Cats, http://www.petmap.com/study_abstracts.html (accessed 10 November 2019).
43. Waddell LS, Brown AJ. Hemodynamic monitoring. In Silverstein DC, Hopper K (eds): *Small Animal Critical Care Medicine*, 2nd ed. St. Louis: Elsevier, 2015, pp 957-962.
44. Cooper E. Hypotension. In Silverstein DC, Hopper K (eds): *Small Animal Critical Care Medicine*, 2nd ed. St. Louis: Elsevier, 2015, pp 46-50.
45. da Cunha AF, Saile K, Beaufriere H, et al. Measuring level of agreement between values obtained by directly measured blood pressure and ultrasonic Doppler flow detector in cats. *J Vet Emerg Crit Car* 2014;24(3):272–8.
46. Jenkins TL, Coleman AE, Schmiedt CW, et al. Attenuation of the pressor response to exogenous angiotensin by angiotensin receptor blockers and benazepril hydrochloride in clinically normal cats. *Am J Vet Res* 2015; 76: 807-813.
47. Lin CH, Yan CJ, Lien YH, et al. Systolic blood pressure of clinically normal and conscious cats determined by an indirect Doppler method in a clinical setting. *J Vet Med Sci* 2006; 68: 827-832.
48. Dos Reis GF, Nogueira RB, Silva AC, et al. Spectral analysis of femoral artery blood flow waveforms of conscious domestic cats. *J Feline Med Surg* 2014; 16: 972-978.
49. Brown SA. Hypertension. International Renal Interest Society. <http://iris-kidney.com/education/hypertension.html> (2016, accessed 14 October 2019).
50. Taylor SS, Sparkes AH, Briscoe K, et al. ISFM Consensus Guidelines on the Diagnosis and Management of Hypertension in Cats. *J Feline Med Surg* 2017; 19(3), 288–303.
51. Acierno MJ, Brown S, Coleman AE, et al. ACVIM consensus statement: Guidelines for the identification, evaluation, and management of systemic hypertension in dogs and cats. *J Vet Intern Med* 2018; 32: 1803– 1822.
52. Chetboul V, Reynolds BS, Trehieu-Sechi E, et al. Cardiovascular effects of dietary salt intake in aged healthy cats: a 2-year prospective randomized, blinded, and controlled study. *PLoS One* 2014; 9: e97862.

Study	Method	n	Systolic	Diastolic	Mean	Machine
<i>Intra-arterial:</i>						
Belew <i>et al.</i> (1999) ²⁷	conscious	13	126 ± 9	106 ± 10	91 ± 11	Dinamap 8300 (Critikon, USA)
Mishina <i>et al.</i> (2006) ²⁸	conscious	20	118 ± 11	78 ± 9	95 ± 10	Model TA11 PA-C40 (Data Sciences International, USA)
Slingerland <i>et al.</i> (2008) ²⁹	conscious	21	132 ± 9	115 ± 8	96 ± 8	Gabarith (Becton Dickinson BV, NL)
Jenkins <i>et al.</i> (2014) ⁴⁶	anaesthetised	6	160 ± 12	138 ± 11	116 ± 8	Model TA11 PA-C40 (Data Sciences International, USA)
<i>Oscillometry:</i>						
Bodey <i>et al.</i> (1998) ³⁵	conscious	104	139 ± 27	99 ± 27	77 ± 25	Dinamap 1846SX (Critikon, USA)
Mishina <i>et al.</i> (1998) ³⁴	conscious	60	115 ± 10	96 ± 12	74 ± 11	USM-700GTM (Ueda Electronic Works, Japan)
Haberman <i>et al.</i> (2004) ²⁰	conscious & anaesthetised	13	132 ± 24	111 ± 24	96 ± 22	Dinamap 8300 (Critikon, USA)
<i>Doppler ultrasonography:</i>						
Sparkes <i>et al.</i> (1999) ¹⁹	conscious	50	162 ± 19			Parks model 811-BTS (Parks Medical, USA)
Haberman <i>et al.</i> (2004) ²⁰	conscious & anaesthetised	13	161 ± 46			Parks model 811 (Parks Medical, USA)
Lin <i>et al.</i> (2006) ⁴⁷	conscious	53	134 ± 16			Parks model 811-B (Parks Medical, USA)
Dos Reis <i>et al.</i> (2014) ⁴⁸	conscious	30	125 ± 16			DV2000 (Medpej, Brazil)
Chetboul <i>et al.</i> (2014) ⁵²	conscious	20	151 ± 5			Parks model 811-BL (Parks Medical, USA)
Payne <i>et al.</i> (2017) ¹⁸	conscious	780	121 ± 16 (IQR 110-132)			Parks model 811-B (Parks Medical, USA)

512 **Table 1:** Published arterial blood pressure (mmHg) values obtained from cats.

	Signalment	Condition	Cuff size	Cuff location	Doppler (mmHg)	PetMAP+ mode A (mmHg)	PetMAP+ mode B (mmHg)
1	17Y BEN FN	Healthy	2.0	L Forearm	170*	209/160 (178)	183/134 (149)
2	11Y BEN M	Vomiting	3.0	L Forearm	160*	156/104 (124)	182/126 (146)
3	8Y DSH MN	Healthy	3.0	L Forearm	139	138/92 (109) *	130/98 (108)
4	3Y DSH FN	Healthy	2.5	R Forearm	146*	181/124†	189/99†
5	10Y DSH MN	HCM	2.5	L Forearm	152	171/113 (132) *	163/114 (130)
6	3Y DSH FN	Lethargy	3.0	L Forearm	131*	159/108 (126)	148/101 (118)
7	10Y DSH MN	Jaundice	2.5	L Forearm	120	157/107 (128) *	159/99 (118)
8	8Y DSH MN	HCM	3.0	R Forearm	148*	150/118 (133)	160/104 (124)
9	13Y DSH FN	Skin wound	2.5	L Forearm	141	189/127 (149) *	175/131 (144)
10	5Y DSH MN	IMHA	2.5	L Forearm	148*	167/82 (107)	153/98 (116)
11	4Y DSH FN	Healthy	2.5	L Forearm	156	199/129 (155) *	185/99 (129)

12	3Y DSH M	RTA	3.0	R Forearm	139*	163/116 (122)	159/104 (120)
13	12Y DSH MN	Anaemia	2.5	R Forearm	147	142/104 (118) *	151/115 (124)
14	4Y DSH FN	RTA	3.0	L Forearm	130*	123/53 (88)	113/65 (84)
15	Adult DSH	Healthy	3.0	R Forearm	132*	139/85 (120)	155/76 (117)
16	9M BEN F	Pelvic fracture	2.5	L Forearm	147	175/127 (144) *	166/116 (134)
17	11M DSH M	Patellar lux.	3.0	R Forearm	141*	181/114 (139)	131/85 (102)
18	6Y DSH MN	Vomiting	2.5	L Forearm	135	156/75 (110) *	162/113 (129)


Table 2: Average of BP measurements obtained from the 18 cats included in the study. (*Y* - years; *M* - months; *BEN* – Bengal; *DSH* – Domestic Shorthair; *F* - female; *FN* – female neutered; *M* - male; *MN* – male neutered; *HCM* – hypertrophic cardiomyopathy; *IMHA* – immune-mediated anaemia; *RTA* – road traffic accident; *L* - left; *R* - right; *Mode A* - Proprietary Optimisations (PPO); *Mode B* - Optimized None (ON); * - used first; † - MAP not recorded)

Device	Blood Pressure	Mean \pm standard deviation (SD)	Bias (mmHg)	Precision (mmHg)	Limits of Agreement (mmHg)	Range* (mmHg)
PetMAP+ mode A	SAP	164.17 \pm 22.36	+21.28	17.67	-13.36 to 55.92	72 to 243
PetMAP+ mode B	SAP	159.17 \pm 20.25	+15.72	16.61	-16.83 to 48.28	50 to 245
Doppler	SAP	143.44 \pm 11.99				110 to 190
PetMAP+ mode A	MAP	128.35 \pm 20.96	-14.94	16.96	-48.18 to 18.30	55 to 195
PetMAP+ mode B	MAP	123.06 \pm 16.24	-20.24	12.72	-45.17 to 4.70	30 to 176

Table 3: A summary of Bland – Altman analyses of PetMAP+ Proprietary Optimisations (PPO) (mode A) and Optimized None (ON) (mode B) function as compared to the Doppler in conscious cats. * - the lowest and highest measurement recorded.

Figure 1:

Oscillometric mode A (SAP) vs Doppler




527

528 Bland-Altman plot of the differences between systolic arterial blood pressures (SAP) measured
 529 by the Doppler and PetMAP+ Proprietary Optimisations (PPO; Mode A) oscillometric device;
 530 measurements from the two devices were plotted against the mean of the Doppler SAP
 531 measurements. Dashed line indicates the bias (average difference), dotted lines indicate upper
 532 limit of agreement (ULA) and lower limit of agreement (LLA).

533

534

535 **Figure 2:**




536

537 Bland-Altman plot of differences between systolic arterial blood pressures (SAP) measured by
538 the Doppler and PetMAP+ Optimized None (ON; Mode B) oscillometric device; measurements
539 from the two devices were plotted against the mean of the Doppler SAP measurements. Dashed
540 line indicates the bias (average difference), dotted lines indicate upper limit of agreement
541 (ULA) and lower limit of agreement (LLA).

542

543 **Figure 3:**




544

545 Bland-Altman plot of the differences between systolic arterial blood pressures (SAP) measured
546 by the Doppler and mean arterial pressure (MAP) by PetMAP+ Proprietary Optimisations
547 (PPO; Mode A) oscillometric device; measurements from the two devices were plotted against
548 the mean of the Doppler SAP measurements. Dashed line indicates the bias (average
549 difference), dotted lines indicate upper limit of agreement (ULA) and lower limit of agreement
550 (LLA).

551

552

553 **Figure 4:**



554 Bland-Altman plot of differences between systolic arterial blood pressures (SAP) measured by
555 the Doppler and mean arterial pressure (MAP) by PetMAP+ Optimized None (ON; Mode B)
556 oscillometric device; measurements from the two devices were plotted against the mean of the
557 Doppler SAP measurements. Dashed line indicates the bias (average difference), dotted lines
558 indicate upper limit of agreement (ULA) and lower limit of agreement (LLA).
559
560